

報告編號： TF15020120

試驗報告

土城場區

發行日期： 2026/03/03

收件/試驗起始日期： 2026/02/12

報告頁次/頁數： 1/3

報告抬頭：大謙科技材料股份有限公司(U3795)

地址：苗栗縣苑裡鎮社苓里社苓 36 之 9 號

試件類別：口罩

數量：1 件

來文字號：--

依委託者所提供來樣資料為：

樣品 1：大謙醫療口罩(未滅菌)



試驗結果：

試驗項目	試驗結果	試驗方法
	樣品 1	
游離甲醛含量 (mg/kg)	N.D.	CNS 14774 T5017- 2022 9.7/CNS 15580-1 L3271-1-2012

註：

1. N.D.=Not Detected.(<16 mg/kg)/未檢出(低於偵測極限值)。
2. 標準限值： ≤ 75 mg/kg。

試驗項目	試驗結果	試驗方法
	樣品 1	
可遷移性螢光物質試驗	無	CNS 14774 T5017- 2022 9.9/CNS 11820 P3082-2007

附記：1.本報告僅對樣品負責，樣品保留期限為一個月。
2.本報告非經本檢測及驗證部書面同意，不得摘錄複製，但全部複製除外。
3.試驗報告所載事項不得作為公開廣告及商業推銷之用。

財團法人紡織產業綜合研究所檢測及驗證部 新北市土城區承天路6號
Department of Testing and Certification, Taiwan Textile Research Institute
No.6, Chengting Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R. O. C.)
Tel: +886-2-22670321 ext. 7107, 7110
Fax: +886-2-22675108, +886-2-22689839

財團法人紡織產業綜合研究所
所長授權核發人

檢測及驗證部 主任 高瑞宏



TF15020120

試驗報告

土城場區

報告編號： TF15020120

發行日期： 2026/03/03

收件/試驗起始日期： 2026/02/12

報告頁次/頁數： 2/3

試驗項目		試驗結果		試驗方法
		樣品 1		
		方法 T	方法 D	
*偶氮染料 (mg/kg)	對-胺基聯苯	N.D.	N.D.	CNS 14774 T5017-2022
	聯苯胺	N.D.	N.D.	9.8/CNS 16113-1 L3280-
	對-氯-鄰-甲苯胺	N.D.	N.D.	1:2019
	2-奈胺	N.D.	N.D.	9.1(方法 D) · 9.2(方法
	鄰-胺基偶氮甲苯	N.D.	N.D.	T) GC-MS
	5-硝基-鄰-甲苯胺	N.D.	N.D.	
	4-氯苯胺	N.D.	N.D.	
	2,4-二氨基苯甲醚	N.D.	N.D.	
	4,4'-二胺基二苯甲烷	N.D.	N.D.	
	二氯聯苯胺	N.D.	N.D.	
	二甲氧基聯苯胺	N.D.	N.D.	
	鄰-二甲基聯苯胺	N.D.	N.D.	
	3,3'-二甲基-4,4'-二氨基 二苯甲烷	N.D.	N.D.	
	2-甲氧基-5-甲基苯胺	N.D.	N.D.	
	4,4'-亞甲雙(2-氯苯胺)	N.D.	N.D.	
	4,4'-二胺基二苯醚	N.D.	N.D.	
	4,4'-硫二苯胺	N.D.	N.D.	
	鄰-甲苯胺	N.D.	N.D.	
	2,4-甲苯二胺	N.D.	N.D.	
	2,4,5-三甲基苯胺	N.D.	N.D.	
鄰-甲氧基苯胺	N.D.	N.D.		
4-胺基偶氮苯	N.D.	N.D.		
2,4-二甲苯胺	N.D.	N.D.		
2,6-二甲苯胺	N.D.	N.D.		

註: 1. "*"委外測試。 2. N.D.=Not Detected.(<5 mg/kg)/未檢出(低於偵測極限值)。 3. 標準限值 ≤ 30 mg/kg。





Report No.: TF15020120

TEST REPORT

TUCHENG

Date of Issue: 03 Mar 2026

Date of Receipt and Test Start: 12 Feb 2026

Page Order/Pages: 1/3

Report Title : Dacian Technology Material Co., Ltd.(U3795)

Address : No. 36-9, Sheling, Yuanli Township, Miaoli County 358 Taiwan

Item : Mask

Quantity : 1 PC

Ref. No. : --

Sample description is given by the client :

Sample 1 : Dacian Medical Face Mask (Non-Sterile)



Test Results :

Test Items	Test Results	Test Method
	Sample1	
Free Formaldehyde Content (mg/kg)	N.D.	CNS 14774 T5017-2022 9.7/CNS 15580-1 L3271-1-2012

Note: N.D.=Not Detected.(<16 mg/kg).

Note: Limit value: ≤ 75 mg/kg.

Test Items	Test Results	Test Method
	Sample1	
Test for Migratable Fluorescent Substances	No	CNS 14774 T5017-2022 9.9/CNS 11820 P3082-2007

Authorized by-president of Taiwan Textile Research Institute
Director, Department of Testing and Certification

Jui-hung kao

Note:

1. This report is only responsible for the submitted sample(s), which will be kept for one month period.
2. This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.
3. The test report should not be used for public advertisement and commercial promotion.

Department of Testing and Certification, Taiwan Textile Research Institute
No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R. O. C.)
Tel: +886-2-22670321 ext. 7107, 7110

Fax: +886-2-22675108, +886-2-22689839



TF15020120

Report No.: TF15020120

TEST REPORT**TUCHENG**

Date of Issue: 03 Mar 2026

Date of Receipt and Test Start: 12 Feb 2026

Page Order/Pages: 2/3

Test Items	Test Results		Test Method	
	Sample 1			
	Method T	Method D		
*AZO DYES (mg/kg)	4-AMINODIPHENYL	N.D.	N.D.	CNS 14774 T5017-
	BENZIDINE	N.D.	N.D.	2022 9.8/CNS 16113-
	4-CHLORO-O-TOLUIDINE	N.D.	N.D.	1 L3280-1:2019
	2-NAPHTHYLAMINE	N.D.	N.D.	9.1(Method D) ·
	6-AMINOAZOTOLUENE	N.D.	N.D.	9.2(Method T) GC-
	2-AMINO-4-NITROTOLUENE	N.D.	N.D.	MS
	4-CHLOROANILINE	N.D.	N.D.	
	2,4-DIAMINOANISOLE	N.D.	N.D.	
	4,4'-DIAMINODIPHENYLMETHANE	N.D.	N.D.	
	3,3'-DICHLOROBENZIDINE	N.D.	N.D.	
	3,3'-DIMETHOXYBENZIDINE	N.D.	N.D.	
	3,3'-DIMETHYLBENZIDINE	N.D.	N.D.	
	3,3'-DIMETHYL-4,4'-DIAMINO-DIPHENYLMETHANE	N.D.	N.D.	
	P-CRESIDINE	N.D.	N.D.	
	4,4'-METHYLENE-BIS-(2-CHLOROANILINE)	N.D.	N.D.	
	4,4'-OXYDIANILINE	N.D.	N.D.	
	4,4'-THIIDIANILINE	N.D.	N.D.	
	O-TOLUIDINE	N.D.	N.D.	
	2,4-TOLUYLENEDIAMINE	N.D.	N.D.	
	2,4,5-TRIMETHYLANILINE	N.D.	N.D.	
O-ANISIDINE	N.D.	N.D.		
P-AMINOAZOBENZENE	N.D.	N.D.		
2,4-Xylidine	N.D.	N.D.		
2,6-Xylidine	N.D.	N.D.		

Note: 1. "*"Subcontracted test. 2. N.D.=Not Detected.(<5 mg/kg) 3. Limit of Standard ≤ 30 mg/kg.

